

The advanced solution for Transcranial Magnetic Stimulation

MagPro X100




- For demanding repetitive protocols and high number of pulses
- Monophasic, biphasic, and half-sine waveforms
- Twin, dual & power theta burst pulses
- Enhanced stimulation capabilities with MagOption
- CE approved and FDA cleared for depression treatment*



Need more options?

- Data storage and reporting capability
- Real-time information on important performance parameters such as remaining stimuli and actual power delivered
- Complete packages available for interleaved TMS/fMRI
- Real double blinded research studies option



Technical data

Mechanical Data		Technical specifications		Ordering Number	
HxWxD, MagPro	210 x 530 x 400mm 8.3 x 20.9 x 15.7 in	Maximum Repetition Rate	100 pps	MagPro X100	9016E0711
HxWxD, MagOption	130 x 530 x 400mm 5.1 x 20.9 x 15.7 in	Pulse mode, X100	Standard	MagPro X100 incl. MagOption	9016E0731
Weight, MagPro	36 kg (79 lbs)	Pulse mode, X100 incl. MagOption	Standard, Dual, Twin and Power		
Weight, MagOption	28 kg (62 lbs)	Waveform, MagPro X100	Monophasic, Biphasic and Biphasic Burst		
		Waveform, MagPro X100 incl. MagOption	Monophasic, Biphasic, Biphasic Burst and Halfsine		
		Current direction	Normal or Reverse		

*Indication for use

MagVenture TMS Therapy® is FDA cleared for:

“Treatment of Major Depressive Disorder in adult patients who have failed to receive satisfactory improvement from prior antidepressant medication in the current episode.”

MagVenture TMS Therapy® is CE approved for:

“Treatment of Major Depressive Disorder in adult patients who have failed to achieve satisfactory improvement from two prior antidepressant medications, at or above the minimal effective dose and duration in the current episode.”

More info at

www.magventure.com

Distributore per l'Italia:



Via Issiglio 95/10, 10141 Torino

+39 011 5821948

+39 011 0433281

info@geasoluzioni.it



Research. Treatment. Results.

©MagVenture A/S • 511-0161, rev. 1.2 • www.magventure.com